DENTOCEM

RESIN CEMENT

DentoCem Resin Cement

F

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DentoCem Resin Cement

GB

Dual cured permanent adhesive resin cement.
Automix.

INDICATIONS:

• Final cementation of crowns and bridges, inlays and onlays, posts and cores, ceramic crowns and Maryland bridges. • Bonding of metal, plastic or ceramic orthodontic attachments. • Adhesive amalgam restoration. • Veneering of alloys. • Implant prosthesis.

PROPERTIES:

• DentoCem Resin Cement is a multi-purpose, dual cured, radiopaque, permanent cement which consists of two easy to mix paste components. • It forms an excellent bond with base alloys, porcelain, dentin, enamel and with amalgams. • The "Automix" system provides a consistent and homogeneous mixture, allows a saving in time and facilitates the use of the product. • DentoCem Resin Cement is a self cured cement which you can light cure for an immediate set and extra strength at the margins. • Light curing DentoCem Resin Cement gives your restorations immediate stability, plus easy cleanup of excess material. • With a simple 10-seconds light precure, you can remove excess cement before the final set. • DentoCem Resin Cement features superior retention and total margin integrity. • When bonding orthodontic brackets, etching of enamel and porcelain is not necessary. • Upon removal of the bracket, no damage is caused to the enamel and porcelain surfaces.

INSTRUCTIONS FOR USE:

Depending on the clinical requirement use the corresponding illustrated card.

Set time: • Working time at room temperature is at least 2 minutes. • Final set time at room temperature is approximately 5 minutes. • Oral set time is approximately 3.5 minutes. • Light cure for 20 seconds.

CAUTION!

DentoCem Resin Cement base and catalyst contain polymerizable monomers which may cause allergic contact dermatitis in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs or if known allergy to methacrylate resin exists, discontinue use.

COMPOSITION:

Base: • Aliphatic polyester urethane acrylate • Hydroxyethyl methacrylate • Crosslinking monomer • Photoinitiators • Adhesive accelerator • Glass filler.

Catalyst: • Triethylenglycol dimethacrylate • Bis-phenol-A-glycidildimethacrylate • Benzol peroxide • Glass filler.

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