# DENTOCROWN

TEMPORARY CROWNS & BRIDGES

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# DENTOCROWN

TEMPORARY CROWNS & BRIDGES

## DentoCrown Temporary Crowns & Bridges

Self-polymerising, paste-paste composite for making temporary crowns and bridges.

### INDICATIONS:

Dentocrown is used to make temporary crowns, bridges, inlays and onlays.

### PROPERTIES:

• High compression resistance: 205 - 215 MPa. • High flexural strength: flexural modulus 1500 - 2500 MPa. / Flexural strength (for 10 min.): 40 - 55 MPa. • Great tensile and fracture strength. Diametrical tension: 25 - 30 MPa. • Easy to remove because of its elastic phase. • Great colour stability. • Easy to finish and polish. • Slightly exothermic polymerisation reaction. • Safety cartridge.


### INSTRUCTIONS FOR USE:

• Insert the Dentocrown cartridge into the mixing gun using a suitable piston. • Remove the cartridge cap and purge the cartridge until you get the same flow rate from both holes. The insert the type-6 mixing nozzle and tighten by twisting clockwise through 90°. The composite is extruded automatically in the correct proportions by pressing on the trigger and can be applied directly. • Do not remove the mixing nozzle from the cartridge after use. Replace the used nozzle with a new one just before re-using. When changing the nozzle, check that the holes are free to ensure that the composite can flow freely.

### APPLICATION:

• Temporary repairs can be made either directly, or with in situ moulds or by using previously taken impressions. • The material for temporary crowns and bridges reaches its elastic phase 1 min to 90 seconds after the start of mixing, i.e. about 30 seconds to 1 min after installing in the mouth. The impression must be removed from the mouth during this phase.

• For the self-etching technique, at a room temperature of 23°C and 50% relative humidity, the working times are as follows:

After about 2-3 minutes remove the material from the impression and dry it (if necessary with a solvent such as ethyl alcohol). This removes the thin inhibition layer caused by oxygen and makes finishing easier. The surface of the temporary repairs must always have a matt, silky appearance. • Carefully remove any excess or rough material in areas close to the gums.

• Finishing and polishing can only be started after the material has completely set (about 4-5 minutes after mixing). • Glue the temporary repair using a eugenol-free, temporary cement (e.g. DentoTemp "ITENA").

### PRECAUTIONS:

• The material for temporary crowns and bridges contains: methacrylates, amines, terebenes, benzoyl peroxide and BHT and must not be used in cases of allergies to one of these compounds. • Before making the temporary repair, the stumps can be isolated to prevent any irritation using fluorinated varnish, stump varnish or similar. • To prevent breakages or injuries caused by fragments, the temporary repair must not be subject to excessive force, especially in the molar region. • Stumps reconstituted using resin compounds (composites) must be suitably isolated so that they don't come apart when the temporary repair is removed (e.g. with glycerine). • Defective hygiene or intensive use of some mouth washes can lead to discolouration.

### STORAGE:

Store at a temperature between 4°C and 23°C. Protect from the light.

# DENTOCROWN

TEMPORARY CROWNS & BRIDGES

## DentoCrown Temporary Crowns & Bridges

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