DENTOTEMP

TEMPORARY CEMENT

DentoTemp Temporary Cement

F

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DentoTemp Temporary Cement

GB

Non-eugenol, acrylic-urethane polymer based temporary cement. Formulated specifically for implants restorations.

INDICATIONS:

- Reliable and quick temporary cementation of temporary crowns and bridges.
- Permanent cementation of implant-retained crowns and bridges for retrievability.

PROPERTIES:

• DentoTemp is unique because it is different than the other temporary or provisional cements on the market. Most current non-eugenol temporary cements are compounds that contain zinc oxide and oil. These materials, while familiar still have many limitations such as clean-up difficulty, interference with permanent cement polymerization and low initial strength. • DentoTemp does not contain essential oils or zinc oxide so it will not harm temporary crowns or interfere with permanent cementation. DentoTemp unique 2-stage polymerization allows for fast and easy clean-up without compromising working time of performance. • DentoTemp is fully compatible with both acrylic and composite (bis-acryl) temporary materials and is indicated on both natural and metal implants abutments. • In addition to being versatile temporary cement, DentoTemp is ideal for permanent cementation of implant-retained crowns. DentoTemp ensures easy cementation and high retention while maintaining retrievability, when desired. • Low solubility in oral fluids and rigid-set mean high retention rate. • Virtually no taste or odor for increased patient comfort.

INSTRUCTIONS FOR USE:

Mixing of DentoTemp:

- 1) Dispense equal amounts of Base and Catalyst onto a mixing sheet.
- 2) Mix together thoroughly with a mixing spatula for 20 to 30 seconds.

Cementation of Temporary Crowns:

Before cementing a temporary crown, check fit, occlusion, proximal contacts and completes all finishing procedures. 1) Dry and isolate the prepared abutments. Absolute dryness is not necessary, but obvious contamination should be avoided. 2) Ensure the interior of the temporary crown is clean and dry. Apply the mixed cement directly into the temporary crown. You have 45-60 seconds working time. 3) Immediately seat the crown and either hold in place or have the patient gently close on a cotton roll until the initial set which is 90-120 seconds. The initial set is characterized by a distinct rubbery phase in which excess cement can be easily removed using an explorer or periodontal probe. 4) Remove excess cement. Take care to pass a probe through interproximal embrasures to eliminate excess cement in those areas. In most cases, excess cement will be extremely easy to remove, usually in several large, easy-to-handle pieces. 5) Final cement set is between 4-5 minutes, depending on conditions.

Note: If cementing a temporary crown over a freshly placed resin core, we recommend very lightly lubricating the core material with a very thin film of a water soluble lubricant to prevent sticking. This step is not absolutely indicated, but certain brands of composite core material may have unpolymerized regions for several hours after curing which may cause sticking.

Cementation of Permanent Implant-Retained Crowns:

The cementation procedure for permanent metal/ceramo-metal crowns over an implant abutment is identical to the cementation of a temporary crown. We do not recommend using a titanium etching agent on the abutments prior to cementation. If removal of the crown is necessary, a gentle tap with a reverse hammer-type crown remover will usually be sufficient.

ESSENTIAL INGREDIENTS:

- Bis-GMA. • Multifunctional methacrylates. • Urethane dimethacrylate. • Polymerization activator.

STORAGE:

- Keep refrigerated. • For dental use only.

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