QUICK BOND

Quick Bond

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QUICK BOND

Quick Bond

Universal, no-rinse, self-priming etchant adhesive system.

GB

INDICATIONS:

• Direct composite restorations. • Adhesive cementation of crow's & bridges, including inlays and onlays.

PROPERTIES:

• Quick Bond is a sixth generation self-etching bonding system, which delivers high bonding values to both, dentin and enamel. • Quick Bond consists of an acidic water base self-etching Primer and a light-cured Adhesive. • It is used to bond composite resin filling material to tooth structure, any self-cure or visible light-cured Composite, Compomer restorative material and with Composite Resin Cements. • Quick Bond is utilizing a moderately acidic Primer. Unlike total-etch, total-bond techniques, the self-etching primer does not open the dentinal tubules completely. The smear layer is solubilized, and due to the high hydrophilic properties of the Primer, it is then able to penetrate into the tubules and peritubular dentin, forming resin tags. In the case of enamel, the Primer produces a significant pattern with enhanced surface area, leading to improved enamel bonding.

INSTRUCTIONS FOR USE:

Application of Quick Bond "A" Prime

1) Apply 3 successive layers of Quick Bond "A" Prime on the prepared, air-dried or blotted tooth surface (time of application : 15 seconds). 2) Leave undisturbed for 5 seconds. 3) Dry gradually until the layer is fixed. 4) Then dry vigorously.

Application of Quick Bond "B" Bond

5) Apply 2 layers of Quick Bond "B" Bond. 6) Remove solvents with a light blast of air. 7) Then dry vigorously during 5 seconds. The surface of the cavity must appear brilliant. 8) Cure for 20 seconds with a visible light cure unit.

When light curing is not an option (bonding in root canal), add equal drops of Bond Activator with Quick Bond "B" Bond in a clean mixing well, combine for 2 seconds and proceed with steps 5, 6 and 7. Restore with direct compositresin of choice. For indirect restorations, use Quick Bond in conjunction with DentoCem Resin Cement.

CAUTION!

Quick Bond contains polymerizable monomers which may cause allergic contact dermatitis in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if you know allergy to methacrylate resin exists, discontinue use.

ESSENTIAL INGREDIENTS:

• Polymerization accelerators • 4-Methacryloxyethyltrimellitic acid • Dendritic Polymer • Water • Acetone • Urethrandiméthacrylate • Urethane modified diméthacrylate • Triéthylenglycoldiméthacrylate • 2-Hydroxyéthylméthacrylate • Multifunctional diméthacrylates • Alkyl Quaternary Ammonium Bentonite • Photoinitiators.

STORAGE:

• Do not store above 24°C • Highly flammable • Contains acetone • Keep refrigerated.

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